

# BITZER Automatic sample analyser

**Automatic sample cleaning  
and determination of  
all relevant grain quality  
values.**

**Quality values  
in 2 minutes**

## Additional benefits

- ▲ Complete laboratory automation for a fast analysis with integrated sample cleaner SLN as well as grain analyser (e. g. Perten or Foss).
- ▲ Time saving by factor 3
- ▲ Automatic process of impurity and quality analysis
- ▲ Storing of all values on one database
- ▲ Results printout for the reference sample optional
- ▲ Integration of further analysing instruments, e. g. falling number device
- ▲ Connection to **BITZER** Web Agrar software or external systems

Filling of grain sample

Touch display for central  
controlling and operating

NIR-Analysis: moisture, protein,  
sedimentation value et al.

Sample cleaner

Print-out of quality values  
for reference sample

