

BITZER Impurity machine MLN

Impurity analysis – safe, clean, fast & quiet

Automatic sample cleaning and determination of the
values of clean grain and impurities via
legal for trade scales in 1 minute



For further information see
www.bitzer-waage.de/en
Simply scan the QR code →



NEW


BITZER Impurity machine

Options:

- ▲ Sieve depot for the storage of different sieves



- ▲ Documentation by reference sample label via printer connection



Print-out of quality values for reference sample

Benefits:

- ▲ Entry level model to improve efficiency in the grain laboratory
- ▲ Time saving due to aspiration, sieving and legal for trade weighing in one process
- ▲ Cost saving due to substitution of separate laboratory scales and sample cleaner (**Pfeuffer Type MLN**)
- ▲ Modular system: expandable to BITZER Automatic sample analyser for complete automation of grain quality determination including connection to grain analyser
- ▲ Comfortable dust protection and sound insulation

Operation panel for central controlling and operating

Filling of grain sample

Sample cleaner



DEALERS
WANTED

BITZER Wiegetechnik GmbH
Benzstraße 3
D-31135 Hildesheim

Tel. +49 5121/ 78 28 -0
Fax +49 5121/7828-78

www.bitzer-waage.de/en
E-Mail: info@bitzer-waage.de

